



PRESS RELEASE: For Immediate Release

Data Innovations Releases Instrument Manager v8.09

May 1, 2009 - SOUTH BURLINGTON, VT: Data Innovations announces the release of Instrument Manager™ (IM) version 8.09, which includes a new Specimen Management Workspace Popup Preferences configuration option, and enhanced Bio-Rad QC OnCall/Unity Real Time Integration. Data Innovations provides free upgrades from previous Windows versions for those enrolled in our support and maintenance programs.

The new Specimen Management Workspace Popup Preferences option provides the ability to select whether to use confirmation windows and dialogue boxes, and to define the default option for these dialogue boxes. This configuration flexibility allows users to streamline the keystrokes needed to perform frequently used actions such as releasing specimens. For actions such as Order Rerun/ Reflex and Order New Tests, the configuration allows users to set the default behavior for which tests should be rerun and for which destination connection.

The Quality Control integration for QC OnCall/Unity Real Time was enhanced to support up to nine levels of controls to be configured for each lot number. This enhanced support for processing any 4 levels of a 9 level quality control lot is provided in advance of the anticipated release of Bio-Rad's Unity Real Time 1.0 service pack 11, slated for the Fall 2009.

Many more enhancements have been made throughout IM. A complete list of these enhancements is included in the IM v8.09 release notes which are accessible from our Customer Web Portal (CWP).

For more information about Instrument Manager v8.09 contact our regional Sales & Marketing Departments.

Region	Sales & Marketing Contact Information
Data Innovations North America	802-264-3470 northamerica-sales@datainnovations.com
Data Innovations Europe	32 2 332 24 13 europa-sales@datainnovations.com
Data Innovations Latin America	55-11-38013283 latinamerica-sales@datainnovations.com
Data Innovations Asia	852-9362 6423 asia-sales@datainnovations.com