



**PRESS RELEASE: For Immediate Release**

**Data Innovations Releases Instrument Manager v8.07**

**July 25, 2008 - SOUTH BURLINGTON, VT:** Data Innovations announces the release of Instrument Manager™ (IM) version 8.07, which includes a new Archiving feature, a re-designed rules interface, and enhanced backup functionality. As always, free upgrades from previous Windows versions are included in all our support and maintenance programs.

Archiving enables the long term storage of data in a separate database for historical and record keeping purposes. Archived data is viewable in Instrument Manager, and can optionally be sorted, cataloged, and stored as PDF files for offline long term storage.

The redesigned graphical user interface for the rules engine, now concurrently accessible by multiple users, gives greater flexibility in creating, maintaining, and testing rules. Rules may be created and edited by free-text typing or by drag-and-drop from a toolbox. Other efficiency features such as a library of rules, frequently used items, and the concept of using value lists within rules have been added.

In order to aid laboratories with regulatory compliance, testing of rules has been greatly enhanced. Users may now build test suites of data, save them, run them through the rules, and save the output. In addition, the rules engine now includes separate test and live rule sets with easy transfer of rules form one environment to the other.

Enhanced functionality allows users to backup set database files manually or automatically at a user defined specified day and time. Backups can be done to any network drive or device without any system downtime.

Many more enhancements have been made throughout IM. A complete list of these enhancements is included in the IM v8.07 release notes, which are attached to this message for your convenience.

For more information about Instrument Manager v8.07 contact our regional Sales & Marketing Departments.

Region	Sales & Marketing Contact Information
Data Innovations <b>North America</b>	802-264-3470 <a href="mailto:sales@datainnovations.com">sales@datainnovations.com</a>
Data Innovations <b>Europe</b>	32 2 332 24 13 <a href="mailto:europe-sales@datainnovations.com">europe-sales@datainnovations.com</a>
Data Innovations <b>Latin America</b>	55-11-38013283 <a href="mailto:latinamerica-sales@datainnovations.com">latinamerica-sales@datainnovations.com</a>
Data Innovations <b>Asia</b>	852-9362 6423 <a href="mailto:asia-sales@datainnovations.com">asia-sales@datainnovations.com</a>

###