

Bio-Rad's Unity Real Time™ software and Data Innovations' Instrument Manager™ (IM) automatically process quality control data on a real-time basis.

This Integrated QC Solution... includes a single point of setup with information entered in the Unity Real Time software available in IM. Unity Real Time and IM can be installed on the same or different computers.

The integration between IM and Unity Real Time provides for real-time analysis of quality control data. When coupled with IM's Specimen Management option, users have the power to define rules that can hold or rerun results based on a data rejection return signal from Unity Real Time. Other functionality, such as entering corrective actions into the Action Log are also integrated.



Not Just QC Software... this offering includes a complete suite of Quality Assurance tools, such as Bench and Supervisor Review, Unity™ Interlaboratory Program, InstantQC™ and QCNet™, as well as the optional Westgard Advisor™, the automatic QC rules selection engine.

Fast Facts: Bio-Rad's Unity Real Time™ Software...

- **Westgard Advisor™** - Recommends and automatically applies optimal QC rules with the use of patented technology. Saves time and money associated with the reduction of unnecessary repeats and troubleshooting.
- **Flexible QC** - Easily set rules for automatic evaluation of QC data, gauge performance with various graphic displays, and access statistical reports for review and archiving.
- **Analytical Goals** - Review the effectiveness of statistical rules. Set target goals for total error, imprecision and comparative bias based on biological and clinical principles.
- **Bench and Supervisor Review** - Display all pending QC data with rule evaluation for both bench and supervisor level reviews, and document review for regulatory compliance.

- **Connectivity** - Automatically import QC data, thereby eliminating manual entry and data entry errors.
- **QCNet™** - Access peer reports, continuing education, software training, control product inserts, Material Safety Data Sheets and diagnostic reference links through a secure website.
- **Unity™ Interlaboratory Program**
Monthly Reports - Comprehensive monthly comparison reports with the industry's largest peer groups.
InstantQC™ - On-demand peer comparison for immediate analysis without deadlines. Results are posted to QCNet™ in easy-to-use tables and graphs.



Data Innovations offers the most comprehensive Middleware system in the market to manage laboratory operations, including pre-analytical, analytical, and post-analytical sample processing and non clinical tasks such as equipment maintenance and specimen archiving. **Instrument Manager™ (IM)**, the Gold Standard in Middleware, is available independent of instrument or LIS vendor, laboratory department, and clinical discipline.

DATA INNOVATIONS, INC.®
 Phone: (802) 658-1955
 Fax: (802) 658-2782
 Email: sales@datainnovations.com

DATA INNOVATIONS EUROPE
 Phone: 32 (0)2 332 24 13
 Fax: 32 (0)2 376 43 84
 Email: europe-sales@datainnovations.com

DATA INNOVATIONS ASIA
 Phone: (852) 2398 3182
 Fax: (852) 2398 8667
 Email: asia-sales@datainnovations.com

DATA INNOVATIONS LATIN AMERICA
 Phone: 55-11-3801-3283
 Fax: 55-11-3871-9592
 Email: latinamerica-sales@datainnovations.com



Bio-Rad Laboratories, Inc. (www.bio-rad.com) is a multinational manufacturer and distributor of life science research and clinical diagnostic products. Bio-Rad is headquartered in Hercules, California, and serves more than 85,000 research and industry customers worldwide through its global network of operations. For more information, please call 800-2-BIO-RAD (800-224-6723) or contact your local Bio-Rad representative.