



## Sample processing for better turnaround time

No matter how many “special circumstances” notes are papering your lab walls and instruments, Data Innovations’ middleware eliminates them all.

With Instrument Manager, all the “if/then” processing instructions ... all the exceptions ... and all the physician requests are written into the software, making your workflow smarter and your results less operator-dependent.

Is your lab **IM**powered™ ?



## Turning STAT requests into STAT results

Our Specimen Routing logic traffics samples through the system based on priority and based on test-load factors. This streamlines the testing to deliver better turnaround times – even when you experience changes in your personnel. You gain the power to balance the workload on the equipment and allocate human resources where they best support your total operation.

- Automatically determines most efficient path for sample processing
- Dynamically alters sample path based on instrument(s) availability/capacity, results, add-ons testing such as repeat/rerun/reflex and more
- Works with automation lines or stand-alone
- Optimize routing to manage carry-over sensitive testing