



WITH STRENGTH AND FLEXIBILITY, YOUR LABORATORY IS IMPOWERED™

## IMSolutions – Hematology™

Today's clinical labs are busier and financial constraints tighter than ever before, thus the need to work smarter and make each man-hour count. With Data Innovations' IMSolutions Hematology,

**This turnkey solution consolidates information for concise viewing providing a holistic approach to results review.**

you get a turnkey middleware solution that eliminates many of the mundane tasks of the past because only "true" exceptions are

presented. In addition, lab technicians can view all the consolidated information on one concise screen. At last, a holistic approach for results review.

The human element is clearly the most important component in the clinical laboratory. When you IMpower your lab you can utilize that resource wisely and efficiently. Data Innovations middleware can help you by:

- Reviewing patient results without the need for technologist intervention
- Relieving tech of mundane tasks, allowing time to attend other duties
- Affording the ability to consolidate multiple results on one screen

### How it works for Centra Health Lab:

*Centra Health has been a huge proponent of the power of middleware to further streamline workflow with sophisticated decision management and rules-based intelligence. By using well-written rule sets for hematology, less than 10% of the CBC's are flagged, allowing for concentrated review of the exceptions. Says Glen McIver with Centra Health, "...the objective is for the lab techs to only have to sign in to one system and on one screen be able to see the cell counter's scattergrams, and on that same screen see the cell counts and their manual keyboard differential tally. And, to be able to do that for any and all analyzers, even if in other buildings. It's the sophisticated stuff that only with the middleware (from DI) will we be able to do."*

### Other features of the IMSolutions Hematology module include:

- Turnkey middleware solution with hardware, software, onsite installation, onsite training and support included
- Support for multiple hematology instruments, LIS connectivity and multiple, concurrent users
- A complete hematology rule set for autoverification that includes the International Society for Laboratory Hematology Consensus Guidelines
- Pre-defined, customizable results review screens that support instrument graphics, previous results, comments and integrated cell counter to perform cell morphology and manual differential counting
- Different configurations are optimized to work with the leading hematology instruments

Data Innovations offers a variety of solutions to your most intractable challenges; solutions that will make your lab run smoother and more efficiently.

Click below to learn more about the Instrument Manager software Solutions

Hematology



With Instrument Manager™, your lab is IMpowered™ to meet demands now and in the future.

For more information, call our solution specialists at 802.264.3470 or visit our website: [www.datainnovations.com](http://www.datainnovations.com).

Document ID: 1629 | Version 2.0

