



Click below to learn more about the Instrument Manager software Solutions

WITH STRENGTH AND FLEXIBILITY, YOUR LABORATORY IS IMPOWERED™

Moving Averages Solution

The reputation of your lab depends on the integrity of your results. Yet, with increased workload and a

Instrument Manager provides continuous monitoring of stability, providing early detection and notification for subtle analytical shifts.

diminishing number of qualified personnel your resources are stretched almost to the breaking point. With Moving Averages from Data Innovations, every sample, every test that comes through

your lab becomes another positive measure of your lab's performance. Your physicians and their patients put their trust in you, Moving Averages ensures that trust is always well placed.

Moving Averages provides a tool to continuously monitor stability, providing early detection and notification for subtle analytical shifts not always detected by standard QC methods. Moving Average Desktop™ allows you to monitor any tests, any work areas from anywhere on your network.

Other features of Moving Averages include:

- User-configurable protocols—define which instruments, tests, number of data points, and algorithms to be used in monitoring performance
- Levey-Jennings style charts to display performance of multiple instruments in real time from anywhere on the network via Moving Averages Desktop
- Users can define the actions to occur automatically when different thresholds are crossed, such as starting to hold tests for verification
- User-defined criteria to automatically exclude results from the calculations and manual exclusions with audit logs

With Data Innovations' unique and powerful Notifier application, notifications can alert you immediately when any work area, any instrument, any test being monitored begins to drift or shift. Along with Moving Averages, Notifier is a key component for any comprehensive quality assurance system. You'll be alerted before a small concern becomes a big problem.

Moving Averages



With Instrument Manager™, your lab is IMpowered™ to meet demands now and in the future.

For more information, call our solution specialists at 802.264.3470 or visit our website: www.datainnovations.com.

Document ID: 1166 | Version 2.0

