

Introducing IM 2020 Edition, version 8.17.00

Learn what's new in Instrument Manager

Your lab is evolving. At Data Innovations, our goal is to help you keep pace. That's why we're proud to announce the release of our newest version of Instrument Manager: IM 2020 Edition, version 8.17.00.

IM 2020 Edition, v8.17.00, includes new features labs can use today to better serve the patients who are at the center of everything you do – and advances the ongoing modernization of the Instrument Manager platform.

Your lab has a lot to gain

Use this IM Upgrade Calculator to learn how much.

What's new?



Timezone awareness



Caché database upgrade



LDAP user auto-provisioning



Specimen Management
Workspace updates

You'll see these new capabilities:

- **Timezone Awareness:** Ability to display and enter results across multiple time zones for labs where IM server and lab users and instruments are not located in the same time zone
- **Caché Database Upgrade:** Upgraded database enables labs to deploy Instrument Manager on Microsoft Windows Server 2019 and benefit from associated Microsoft security updates
- **LDAP user auto-provisioning:** It's even easier to use Instrument Manager in conjunction with LDAP/Active Directory. User access in Instrument Manager can be built off existing LDAP profiles, versus having to configure users manually
- **Specimen Management (SM) Workspace:** We're continuing to enhance Specimen Management. Now, you can launch 3rd party applications – web links or PDF downloads – from comments in SM Workspace

Ready to upgrade to IM 2020 edition,
version 8.17.00?

Contact Us