



LAB GPS™, V2.3
SYSTEM SPECIFICATIONS

Because each customer situation is unique, discuss system requirements with Data Innovations® Sales and Technical teams prior to any purchase. For budgetary and planning considerations, Data Innovations outlined two categories of system specifications to provide guidelines for the System types required to best suit the intended use of Lab GPS™ and your production volume.

Workstation Requirements	Server Requirements
<p>Workstation System Specifications¹</p> <ul style="list-style-type: none"> • 25 Mbps download speed or higher • 5 Mbps upload speed or higher • 4 GB or more of RAM • Dual Core Processor 	<p>Installation Requirements</p> <ul style="list-style-type: none"> • Instrument Manager™ version 8.17 or higher • diietldl version 8.15.0009 or higher • Microsoft® .NET Framework 4.7.2 or higher • Administrative privileges to install driver and Windows service
<p>Hardware</p> <p>Monitor with minimum 1920 x 1080 resolution, keyboard, and mouse.</p>	<p>Supported Microsoft® Operating Systems (64-bit)</p> <ul style="list-style-type: none"> • Windows® 10 Professional • Windows® 11 Professional • Windows Server® 2019² • Windows Server® 2022
<p>Supported Browsers</p> <ul style="list-style-type: none"> • Google Chrome® v136 • Microsoft Edge® v136 <p>NOTE Other browsers, such as Mozilla® Firefox®, may be compatible with Lab GPS™, but they have not been tested <i>and their capabilities and/or restrictions are unknown.</i></p>	<p>Network Specifications</p> <ul style="list-style-type: none"> • Network interface card installed (NIC) • TCP/IP networking installed • Static IP address • Ports 1972, 57772, and 443 for communication with InterSystems Caché/IRIS® • Remote Procedure Call (RPC services) must be running for driver configuration • TLS 1.2 protocol or higher

1.Contact Data Innovations’ Sales by calling (888) 299-1750 or emailing northamerica-sales@datainnovations.com to discuss the specific intended use of Lab GPS™ and production volumes prior to purchasing computer equipment.

2.Microsoft .NET 4.7.2 will need to be installed.

Single Sign-On (SSO) Requirements

Lab GPS is based on SAML 2.0 protocols and supports both Service Provider (SP) and Identity Provider (IdP) initiated SAML SSO flows.

- SP-initiated SSO: A Lab GPS user enters the Lab GPS URL in the browser and the user is redirected to the customer IdP for authentication. Upon successful authentication, IdP redirects user back to Lab GPS. This flow is recommended and enabled by default.
- IdP-initiated SSO: A Lab GPS user first logs in to the IdP and is presented with all the applications the user has access to. The user then selects Lab GPS and is redirected to Lab GPS.

Mobile App Requirements

- Network specifications:
 - TLS 1.2 protocol
- Supported operating systems:
 - Android v15
 - iOS v18

NOTE Other operating system versions that support TLS 1.2 may be compatible with the Lab GPS™ mobile app, but they have not been tested *and their capabilities and/or restrictions are unknown.*
- Single sign-on enabled



Data Innovations LLC

NORTH AMERICA
463 Mountain View Drive, Suite 305
Colchester, VT USA 05446
Fax: (802) 658-2782

Support

Phone: (802) 658-1955
Email: northamerica-support@datainnovations.com

Sales

Phone: (802) 264-3470
Email: northamerica-sales@datainnovations.com

D.I. Hong Kong LTD.

ASIA
Fax: +852 2398 8667

Support

Phone: +852 2398 3182
Email: asia-support@datainnovations.com

Sales

Phone: +852 2398 3182
Email: asia-sales@datainnovations.com

Data Innovations Latin America Ltda.

LATIN AMERICA
Fax: +55 (11) 3871-9592

Support

Phone: +55 (11) 3801-3283
Email: latinamerica-support@datainnovations.com

Sales

Phone: +55 (11) 3801-3283
Email: latinamerica-sales@datainnovations.com

Data Innovations S.A.

EUROPE
Fax: +32 2 376 43 84

Support

Phone: +32 2 332 24 13
Email: europa-support@datainnovations.com

Sales

Phone: +32 2 770 62 22
Email: europa-sales@datainnovations.com

About Data Innovations

Data Innovations is a global software company that is passionate about excellence in patient care. Through innovative solutions and world-class service, we enable enterprise management of hospital and independent laboratories. Founded in 1989 with headquarters in Vermont, Data Innovations serves more than 6,000 hospitals and laboratories in 80+ countries.

